
Acland® Micro
Vascular Clamps
*Acland® Mikro-
Gefäßklemmen*

Acland® Micro Vascular Clamps 7.01-12
Acland® Mikro-Gefäßklemmen

Basic Informations <i>Grund-Information</i>	7.02-06
Acland® Micro Vascular Clamps (single- and double-clamps) <i>Acland® Mikro-Gefäßklemmen (Einzel- und Doppelklemmen)</i>	7.07
Approximator clamps with frame Clamps by Acland®, angulated Single Use Micro Vascular Clamps <i>Approximatorklemme mit Rahmen Klemmen nach Acland®, abgewinkelt Einweg-Mikro-Gefäßklemmen</i>	7.08-09
Clamp Applying Forceps <i>Klemmenanlegepinzette</i>	7.10-11
Clamp boxes <i>Klemmdosen</i>	7.12

S&T[®]
Basic
informations
*Grund-
Information*

Clinical Micro Clamps by Acland[®]

Prof. Dr. Robert D. Acland, University of Louisville,
Louisville. Ky. USA

The Original

What started as a laboratory instrument over thirty years ago, has become the most popular micro clamp due to the many advantages this clamp series has for its users:

- light and compact design
- extremely durable
- corrosion resistant metallurgy
- a unique gripping surface
- smooth sliding bar action on all approximator clamps
- clamp profile with parallel sided jaws
- choice of matte or black finish

Offering a range of 5 sizes, the **S&T**[®] clinical micro clamp series will meet the most exacting needs of small vessel surgery in the diameter range from 0.4 to 5 mm.

Clamp Security

All **S&T**[®] clamps are precisely calibrated to exacting standards. All **S&T**[®] clamps undergo a rigorous QC testing procedure.

S&T[®] clamps are **100% inspected** prior to shipment. **S&T**[®]'s clinical micro clamps are the gentlest and least traumatic of all small vessel clamps.

Important!

Acland[®] is a registered trademark of **S&T**[®] AG.

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Basic informations Grund- Information

Slippage

The flat and non-serrated "Ultra-Grip" surface of these clamps ensures that the pressure is spread evenly over the entire area of the vessel wall. Minimizing over-pressure peaks, which may lead to damage of the vessel interior.

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Options


Clamp size?

In choosing the correct clamp size for your application, please refer to the section "Clamp Size and Clamp Pressure" below.

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"A" or "V" Pattern?

Throughout the series, the "V" pattern is the all-purpose clamp, commonly used on all veins and almost all arteries. The "A" pattern is designed for exceptionally thick-walled and slippery arteries. The slight incurve of the tips of the "A" pattern jaw provides extra security against slippage.

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Black or Plain Finish?

All clamps are now available in either a black or matte (dark gray) finish. The black finish helps minimize eye fatigue, and also helps avoid burned-out highlights in photography and on-screen. The matte finish features a dark gray, minimally reflective surface. The selection of a black or matte finish clamp is a matter of preference.

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Basic informations Grund- Information

Clamp Size and Clamp Pressure

When choosing the right size of clamp for your work, clamp pressure is an important consideration. Clamp pressure is determined by the closing force of the clamp (in grams), divided by the area of the vessel wall that is compressed between the jaws (in square millimeters). Thus, the smaller the area, the higher the pressure is, inversely proportional to the vessel size. Each size of clamp has its own different and carefully calibrated closing force. Each clamp size is designed for a specific range of vessel diameters (Fig. 1).


Fig. 1

Clamp size Klemmengröße	Size range (mm) Gefässdurchmesser (mm)
B-1 / P-1 / ABB-11	0.4 - 1.0 • ○
B-2 / P-2 / ABB-22	0.6 - 1.5 ○ ○
B-3 / P-3 / ABB-33	1.0 - 2.25 ○ ○
HD / HD-D	1.5 - 3.5 ○ ○
RD / RD-D	2.0 - 5.0 ○ ○

Fig. 2

Leakage

Even on the largest, recommended diameter vessel (Fig. 2), each size of clamp will withstand blood pressure within the vessel up to 200 mm Hg, without leakage.

Despite the wide range of closing forces, the pressure exerted by each clamp, big or small, remains consistent throughout the series. Each size of clamp exerts a pressure of 5g/mm² when used on the largest vessel in its range, and 15g/mm² when used on the smallest vessel. Even at the higher range, S&T's clinical micro clamps have the gentlest working pressure of any small vessel clamp available.

The table (Fig. 1) summarizes the relationship between clamp size, recommended vessel diameter and clamp pressure.

"Vessel diameter" refers to the external diameter of the vessel in its natural state of full dilation.

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Basic informations

Grund-Information

There is a 50% overlap in the recommended vessel diameter range, between each clamp size and the next larger or smaller size. Thus, for every vessel diameter, with the exception of those at the extreme ends of the range, there will always be two suitable clamps. Whenever a choice arises, the smaller clamp should be selected to ensure a gentler, less pressured handling of the vessel.

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Our new design

Advantages

- Marking shows the appropriate position of vessels
- Border lines show the range of the vessel diameters
 - The outer lines designate the maximum vessel diameter (Max)
 - The inner lines designate the minimum vessel diameter (Min)
- The **S&T®** logo: The certainty that you are using original Acland® microvascular clamps!

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
Basic
informations
*Grund-
Information*

Single clamps

Single clamps are used in situations where space is at a premium and you find it difficult to place an approximator clamp. Generally you will use single clamps when you are working on end-to-side anastomosis and you need safe haemostasis on the flap vessel.

Packaging unit: 2 single clamps

Approximator clamps

Approximator clamps are the standard clamps for end-to-end anastomosis or if the recipient vessel on an end-to-side anastomosis is to be clamped and approximated. The ease with which vessel ends may be positioned as well as the exceptional access to the procedure site has made these clamps top sellers.

Packaging unit: 1 Approximator clamp

Care of Clamps

Micro clamps are prone to damage from the accumulation of blood clotting in tight places. With proper care and cleaning, this can be avoided. Blood should not be allowed to dry on the clamps during surgery. For cleaning after surgery, we recommend immersion for 30 minutes in hemolytic detergent (please also note the section on "Decontamination"). When immersing the clamps, hold them open for a moment to wet the inner surfaces. Be sure to slide approximator clamps along the bar, back and forth a few times, to clear the sliding units. After immersion, rinse vigorously with a fine, high-pressure water jet to remove any blood residue from spring interspaces, sliding units and suture cleats.

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Acland® Micro Vascular Clamps *Acland® Mikro- Gefäßklemmen*

Approximator clamps with frame

These clamps are generally used in labs or in instances where there is no assistance available. The cleats on the frame are used for belaying the suture material, ensuring perfect positioning during an anastomosis.

Important!

Do not use a clamp on a vessel below the recommended diameter range. Use of any clamp on a vessel too small will result in unacceptably high pressure and risk of vessel damage.

Clamps by Acland®, angulated

These handy clamps are recommended for use in places where straight clamps either are difficult to apply, or obstruct sight in the operative field.

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Acland[®] Micro Vascular Clamps *Acland[®] Mikro- Gefässklemmen*

Applying the Clamps

The three smaller sizes, B-1, B-2, and B-3, are designed to be applied with the Meyer clamp applicator forceps. The forceps are direct acting. When the forceps is squeezed, the clamp opens. When held in the forceps, the clamp can be manipulated firmly in any plane without the danger of slipping or falling out. Clamp applicator forceps are available with or without a lock.

RD and HD clamps are designed to be applied by hand. Their broad, corrugated gripping surfaces make this easy.

Prof. Dr. Viktor E. Meyer, University Hospital Zurich, Zurich, Switzerland

Important!

Applying the three smaller sized B-clamps by hand may damage the calibration. Please use the clamp applying forceps for it.

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Clamp Boxes Klemmdosen

Sterilization case CB-2 for original Acland®-clamps

- Easy access and clearly arranged identification of the clamps
- Suitable for all current sterilization methods
- For instant use with clamp applying forceps
- Facilitates work and classification in OR
- Provides necessary protection during sterilization, transportation and storage
- Economical, durable and functional
- Perforation provides good ventilation to facilitate drying
- A simple, cost-efficient solution to ensure longevity of your original Acland®-clamps!

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