
S&T®

Care of  
Instruments  
*Instandhaltung  
der Instrumente*

0

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# S&T® Care of Instruments *Instandhaltung der Instrumente*

## Instrument Care and Handling

Proper care and handling techniques will help preserve the life of an instrument, and ensure that it functions at peak performance. Avoid having your instruments lay loose together. You should not combine delicate micro instruments with macro instruments as this may cause damage to the delicate tips. Most damage to microsurgical instruments occurs during processing and cleaning, so it is crucial that you handle the instruments with great care. Finally, we recommend you invest in a sturdy instrument container that is specifically designed for the storage, sterilization and transport of microsurgical instruments.

**Important**  
Most instrument damage occurs during processing and cleaning. The best way to avoid this is by handling them with great care. This is especially true for microsurgical instruments.

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# Care of Instruments *Instandhaltung der Instrumente*

## Pre-cleaning

It is best to remove excess tissue or other residues from your instruments as soon as possible, ideally immediately after the procedure. This can be done manually or by an ultrasonic cleaner. We suggest that a disinfecting agent be used. Keep in mind that temperatures below 40°C are needed to avoid protein coagulation.

Make sure to move all approximators, ratchets and sliding parts to ensure that all residues are removed. Clamps, scissors and needleholders need to be opened to allow contact with cleaning solutions.

## Cleaning

In general, you should follow the pre-scribed guidelines established for your institution. It is recommended that cleaning and disinfecting be done in a thermal disinfectant. Initial cleaning commences at temperatures of approximately 40-60°C with pH-neutral or alkaline cleaners. Rinsing is best done with cold water and without additives. Foam generation should be avoided as it may inhibit the cleaning action.

You should follow the manufacturers' instructions regarding use of cleaning agents and cleaning equipment. Water used in all cleaning processes should have a low concentration of chloride ions, as excessive concentrations may lead to pitting and other forms of corrosion.

Final rinsing should be done with distilled / deionized water. This helps to avoid staining of the instruments.

### Important

**To avoid damaging the delicate tips of micro instruments, process them in their own designated tray. Avoid overloading the tray, making sure that micro instruments do not come in contact with macro instruments. Scissors and needleholders should be opened to ensure proper cleaning.**

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# S&T®

## Care of Instruments

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Be mindful of instruments with large surface areas. They can overshadow smaller instruments situated below them, inhibiting them from being cleaned properly and effectively.

Instruments should be dried thoroughly after cleaning, to avoid corrosion.

Only clean instruments should be sterilized.

#### **Important**

**New instruments should be processed and sterilized prior to their first use.**

### **Disinfection**

For the disinfecting process, we recommend using a cycle with the highest temperatures (usually 80-95°C).

### **Inspection**

Instruments should be inspected after each use and immediately after cleaning. The instruments tips should approximate well. Ensure that all scissors and needleholders are in proper working order. Pay special attention to hinges and box locks. Make sure that they are thoroughly cleaned and function smoothly. Check for bent and broken tips. Instruments that are damaged should be removed from service and either be repaired or replaced.

Instruments that are beyond repair should be discarded, ensuring that they will not find their way back into the operating room.

If corrosion is apparent, the instrument should be removed from service and sent out for repair. Depending on the nature of the corrosion, some mild forms may be polished away. In the case of serious pitting or corrosion, the instrument should be discarded as to avoid cross contamination with other instruments.

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# Care of Instruments *Instandhaltung der Instrumente*

## Lubrication

**Instruments with moving parts perform best when properly lubricated.**

Instrument lubrication not only makes an instrument perform better, but also reduces the effects of corrosion, ultimately prolonging the life of the instrument. We recommend any paraffin-based lubricant that is certified for use with medical devices. Remember to use only enough to wet the relevant areas, allowing instruments to drip dry before placing them in a tray or wrapping.

**S&T®** instruments have a **lifetime warranty** on workmanship and materials. Instruments that are judged to be defective should be returned to **S&T®**. Contact your local representative. When returning an instrument, supply as much of the following information as possible: Article code or catalog number; numerical code, located next to the CE mark; and description of the defect. This information is used to improve **S&T®** instruments and ensure that our products remain state-of-the-art.

## Packaging

Packaging plays an important role in the sterility of your instruments. Corrosion may occur if instruments are handled improperly during the packaging process, as the drying cycle in the sterilizer may be inhibited. Further, it is possible to invalidate the entire sterilization process if defective packaging materials do not properly protect the sterilized instruments.

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# Care of Instruments *Instandhaltung der Instrumente*

## Sterilization

In general, you should follow your institution's set guidelines and procedures when sterilizing microsurgical instruments. The most common form of sterilization is steam autoclaving. Microsurgical instruments should be sterilized in their own case. We recommend autoclaving your **S&T®** instruments in saturated steam, at a minimum of 132-134°C, for 3 -18 minutes. Stained instruments are an indication of problems having to do with cleaning/sterilization processes relating to one or more of the following:

- cleaning agents
- steam quality in sterilization
- water quality
- other agents

### Important

**Non-stainless or corroded instruments may damage other instruments during sterilization.**

## Storage

Storage sites should be clean, dust-free and well ventilated. Environmental parameters such as temperature, humidity and radiation should be congruent with the standards for human habitation. Humidity must be below the dew point, at all times.

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## Care of Instruments

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#### Use

The life of an instrument may be extended if it is cleaned as often as possible during a procedure. Blood and other contaminants should be removed from the surface of the instrument before it is allowed to dry. This is best done with distilled /deionized water, **not** with saline or ringers solution.

Instruments that appear or become damaged during a procedure should be removed from use and marked for later evaluation, repair and/or replacement.

An influencing factor in an instrument's functionality is the transportation of the instrument between processing facilities, storage and use. We recommend the user evaluate methods and establish a system that is as simple and safe as possible.

#### Important

The guidelines that we have provided here are general in nature. Care and handling procedures vary from country to country and from institution to institution. You should adhere to procedures established by your local protocol. If you require additional information please contact us directly or your local distributor or visit also the website of the "Instrument Preparation Working Group" on [www.a-k-i.org](http://www.a-k-i.org) for a free pdf-download of the "Red Brochure" (Guidelines for proper maintenance of instruments).

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