


Peanger | Artery forceps | Arterienklemmen |  
Pinze emostatiche | Pinzas hemostáticas | 止血鉗子



Uniform, rounded tips.

Bevelled edges to protect the tissue.

Box joints ensure perfect jaw alignment. No sharp edges.

Ratchets ensure effortless and safe locking and unlocking.

Shanks are slim and flexible for easy handling.


Deeply cut, well defined serrations, mesh perfectly.



### Pean STILLE peang | Pean STILLE Artery Forceps

General utility forceps used in many surgical disciplines.

103-10015-131 Straight 13.5 cm 5.25"



### Halsted-Mosquito STILLE peang | Halsted-Mosquito STILLE Artery Forceps

Hemostat without teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.



103-10030-121



103-10032-121



103-10032-181

103-10030-121 Straight 12 cm 4.75"

103-10032-121 Curved 12 cm 4.75"

103-10032-181 Curved 18 cm 7"



### Halsted-Mosquito STILLE gracil peang | Halsted-Mosquito STILLE delicate Artery Forceps

Delicate hemostat without teeth for small operating field.




103-10026-121



103-10028-121

103-10026-121 Straight 12 cm 4.75"

103-10028-121 Curved 12 cm 4.75"



### Halsted-Mosquito STILLE peang | Halsted-Mosquito STILLE Artery Forceps

Hemostat with 1x2 teeth, used in small operative fields. For plastic surgery, neuro, hand or any delicate surgery.



103-10035-121



103-10037-121

103-10035-121 1x2 teeth Straight 12 cm 4.75"

103-10037-121 1x2 teeth Curved 12 cm 4.75"

### Kocher STILLE peang | Kocher STILLE Artery Forceps

Hemostat with 1x2 teeth. Widely used when a short instrument is needed. Can also be used as tissue forceps to retract or lift an organ.



103-10115-141


103-10115-141	1x2 teeth	Straight	14 cm 5.5"
103-10117-141	1x2 teeth	Curved	14 cm 5.5"



### Kelly STILLE peang | Kelly STILLE Artery Forceps

Fairly short hemostat with jaws, half-serrated.

103-10150-141	Straight	14 cm 5.5"
103-10152-141	Curved	14 cm 5.5"






**Kelly-Rankin STILLE peang | Kelly-Rankin STILLE Artery Forceps**


Longer version of the Kelly hemostat with jaws half-serrated.

103-10160-161	Straight	16 cm 6.25"
103-10162-161	Curved	16 cm 6.25"



**Christophe-STILLE peang | Christophe-STILLE Artery Forceps**


Curved hemostat without teeth.



103-10172	Curved	16 cm 6.25"
-----------	--------	-------------

### Crile STILLE peang | Crile STILLE Artery Forceps

General purpose hemostat without teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.




103-10180-141	Straight	14 cm 5.5"
103-10180-161	Straight	16 cm 6.25"
103-10180-181	Straight	18 cm 7"
103-10182-141	Curved	14 cm 5.5"
103-10182-161	Curved	16 cm 6.25"
103-10182-181	Curved	18 cm 7"



### Crile STILLE peang | Crile STILLE Artery Forceps

General purpose hemostat with 1x2 teeth. Used in relatively small operating fields, when more strength is needed than a Mosquito can offer.




103-10185	1x2 teeth	Straight	14 cm 5.5"
103-10185-161	1x2 teeth	Straight	16 cm 6.25"
103-10187-141	1x2 teeth	Curved	14 cm 5.5"
103-10187-161	1x2 teeth	Curved	16 cm 6.25"



### Leriche STILLE peang | Leriche STILLE Artery Forceps

Fairly tapered hemostat. For use in small operative fields.




103-10190		Straight	15 cm 6"
103-10192		Curved	15 cm 6"
103-10195	1x2 teeth	Straight	15 cm 6"
103-10197	1x2 teeth	Curved	15 cm 6"



### Rochester-Péan STILLE peang | Rochester-Péan STILLE Artery Forceps

Sturdy hemostat. Relatively blunt tips. Used in all types of surgery.




103-10200-161		Straight	16 cm 6.25"
103-10200-181		Straight	18 cm 7"
103-10202-161		Curved	16 cm 6.25"
103-10202-181		Curved	18 cm 7"



### Rochester-Ochsner STILLE peang | Rochester-Ochsner STILLE Artery Forceps

Sturdy hemostat with 1x2 teeth. Used when a strong forceps is needed, both to occlude blood vessels and to lift or retract organs.




103-10205-161	1x2 teeth	Straight	16 cm 6.25"
103-10205-181	1x2 teeth	Straight	18 cm 7"
103-10207-161	1x2 teeth	Curved	16 cm 6.25"
103-10207-181	1x2 teeth	Curved	18 cm 7"



### Rochester-Péan STILLE peang | Rochester-Péan STILLE Artery Forceps

Sturdy hemostats. The 20 cm has short and the 24 cm has long jaws. Also referred to as hysterectomy forceps.




103-10210-201	Straight	20 cm 8"
103-10210-241	Straight	24 cm 9.5"
103-10212-201	Curved	20 cm 8"
103-10212-241	Curved	24 cm 9.5"
103-10212-301	Curved	30 cm 12"



### Rochester-Ochsner STILLE peang | Rochester-Ochsner STILLE Artery Forceps

Sturdy hemostat with 1x2 teeth. Used when a strong and long forceps is needed, both to occlude blood vessels and to lift or retract organs.




103-10215-201	1x2 teeth	Straight	20 cm 8"
103-10215-241	1x2 teeth	Straight	24 cm 9.5"
103-10217-201	1x2 teeth	Curved	20 cm 8"
103-10217-241	1x2 teeth	Curved	24 cm 9.5"



### Bengolea STILLE peang | Bengolea STILLE Artery Forceps

A long hemostat with straight or curved jaws without teeth.



103-10220-201	Straight	20 cm 8"
103-10222-201	Curved	20 cm 8"





### Olivecrona STILLE peang | Olivecrona STILLE Artery Forceps

16 cm straight hemostat. Useful as a long hemostat with fine points.



103-10240-161

103-10240-161	Straight	16 cm 6.25"
---------------	----------	-------------




### Adson STILLE peang | Adson STILLE Artery Forceps

Popular model when a little longer and fairly delicate hemostat is needed.

103-10260-191	Straight	19 cm 7.5"
---------------	----------	------------

103-10262-191	Curved	19 cm 7.5"
---------------	--------	------------



### Dandy STILLE peang | Dandy STILLE Artery Forceps

Side curved hemostat without teeth. Often used in neuro surgery.



103-10310

103-10310 14 cm 5.5"



### Schmidt STILLE peang | Schmidt STILLE Artery Forceps

A hemostat popular in general surgery.




103-10411-201

103-10411-201 Curved 20 cm 8"



### Sawtell STILLE peang | Sawtell STILLE Artery Forceps

A hemostat with half-serrated and fairly strongly curved jaws.




103-10413-211	Curved	21 cm 8.25"
103-10413-241	Curved	24 cm 9.5"



### Heiss STILLE peang | Heiss STILLE Artery Forceps

Dissecting and ligature clamps with half-serrated, jaws.




103-10421-211	Straight	21 cm 8.25"
103-10422-211	Curved	21 cm 8.25"
103-10420-201	Strongly curved	20 cm 8"



### Overholt STILLE peang | Overholt STILLE Dissecting Forceps

Strongly curved dissecting forceps.




103-10425	Strongly curved	14 cm 5.5"
-----------	-----------------	------------




### Crafoord STILLE peang | Crafoord STILLE Artery Forceps

Designed for thoracic and cardiovascular surgery. The fine tips make this a very useful instrument. Offers the surgeon Mosquito type tips, but with a length that allows working in a deep cavity.



103-10416




103-10430-241

103-10416	Curved	18 cm 7"
103-10430-241	Curved	24 cm 9.5"



### STILLE artär- och ligaturklämmare | STILLE Artery and Ligature Clamp

Four very delicate models of dissecting forceps with different jaw curvature.



103-10423-1	Slightly curved	20 cm 8"
103-10423-2	Curved	20 cm 8"
103-10423-3	Strongly curved	19.5 cm 7.75"
103-10423-4	Angled 90°	19 cm 7.5"



### STILLE-Ehrén ligaturklämmare | STILLE-Ehrén Ligature Clamp

For fine and precise handling of ligatures in very narrow areas, such as paediatric surgery. The jaws close only at the tips to minimize the risk for accidental clamping of fragile vessels etc. Developed together with Henrik Ehrén, M.D. at Karolinska University Hospital. Matted finish.

103-10428-151 Delicate Strongly curved 15 cm 6"

