



# Klämmare | Clamps | Klemmen | Pinze | Pinzas | 血管鉗子および組織把持鉗子




### STILLE colonklämmare | STILLE Colon Clamp

Colon clamp with atraumatic DeBakey jaws. Long ratchet for precise adjustment.




107-10565 Atraumatic 5 cm jaws 23.5 cm 9.25"



### Wikström STILLE galltång | Wikström STILLE Gall Duct Clamp

Gall duct clamp with angled, serrated jaws.




107-10571-241

107-10571-241 Angled 24 cm 9.5"



### Uppsala STILLE galltång | Uppsala STILLE Gall Duct Clamp

Gall duct clamp with S-shaped, serrated jaws.



107-10581-241

107-10581-241 S-shaped 24 cm 9.5"



**STILLE mini universalklämmare |  
STILLE mini Multi-purpose clamp**

Multipurpose vascular mini clamps with spoon shaped atraumatic DeBakey jaws.



107-10655



107-10656

107-10655	Small	20 mm jaws	15.5 cm 6.25"
107-10656	Angled shanks	20 mm jaws	15 cm 6"



**STILLE mini universalklämmare |  
STILLE mini Multi-purpose clamp**


Multipurpose vascular mini clamps with atraumatic DeBakey jaws.



107-10665



107-10670




107-10671

107-10665	Straight shanks	25 mm jaws	15.5 cm 6.25"
107-10670	Angled 30°	30 mm jaws	15.5 cm 6.25"
107-10671	Angled 50°	30 mm jaws	15 cm 6"



STILLE universalklämmare | STILLE  
 Multi-purpose clamp

Multi purpose clamps with angled, atraumatic DeBakey jaws.




107-10825	Atraumatic	6 cm jaws	23 cm 9"
107-10826	Atraumatic	7.5 cm jaws	22.5 cm 8.75"



**STILLE universalklammare | STILLE  
Multi-purpose clamp**

Multi purpose clamps with angled, atraumatic DeBakey jaws.




107-10872	Atraumatic	6 cm jaws	18,5 cm 7.25"
107-10873	Atraumatic	3.5 cm jaws	15,5 cm 6.25"
107-10876	Atraumatic	8 cm jaws	20 cm 8"




STILLE mini tangential occlusionsklämmare |  
 STILLE mini Tangential Occlusion Clamp


Tangential occlusion mini clamps with three different jaw sizes.  
 With atraumatic DeBakey jaws.



107-10695



107-10696




107-10697

107-10695	Small	14 mm jaws	14 cm 5.5"
107-10696	Medium	15 mm jaws	15 cm 6"
107-10697	Large	24 mm jaws	15 cm 6"




**Satinsky STILLE tangential occlusionsklämmare**  
 | Satinsky STILLE Tangential Occlusion Clamp


Tangential occlusion clamps with atraumatic DeBakey jaws.



107-10840



107-10841



107-10842


107-10840	Atraumatic	Small	22 cm 8.75"
107-10841	Atraumatic	Medium	22.5 cm 9"
107-10842	Atraumatic	Large	23 cm 9"






## STILLE perifer kärlklämmare | STILLE Peripheral Vascular Clamp

Peripheral vascular clamp with atraumatic DeBakey jaws.




107-10871

107-10871 Atraumatic 6 cm jaws 20 cm 8"




**STILLE anastomosklämmare | STILLE  
Anastomosis Clamp**


Anastomosis clamps with atraumatic DeBakey jaws.



107-10860



107-10861




107-10862



107-10860	Atraumatic	1 cm jaws	17,5 cm 7"
107-10861	Atraumatic	1,5 cm jaws	17,5 cm 7"
107-10862	Atraumatic	2 cm jaws	18,5 cm 7.25"

**Lambert-Kay STILLE aorta anastomosklämmare |  
Lambert-Kay STILLE Aortic-Anastomosis Clamp**

Aortic anastomosis clamp with atraumatic DeBaKey jaws.




107-10845

107-10845 Atraumatic 22 cm 8.75"




**STILLE DeBakey aorta aneurysmklämmare** |  
**STILLE DeBakey Aortic-Aneurysm Clamp**

Aortic aneurysm clamp with atraumatic DeBakey jaws.



107-10846




107-10849

107-10846	Atraumatic	8 cm jaws	26 cm 10.25"
107-10849	Atraumatic	7 cm jaws	28 cm 11"



### Ericsson-STILLE carotisklämmare | Ericsson-STILLE Carotid Clamp

Carotid clamp with angled atraumatic DeBakey jaws.




107-10880

107-10880 Atraumatic 17,5 cm 7"



**Ericsson-STILLE profunda arteria femoralis klämmare** | Ericsson-STILLE Profunda Femoral Artery Clamps

Profunda femoral artery clamp with atraumatic DeBakey jaws.



107-10898

107-10898 Atraumatic 18 cm 7"



### Backhaus STILLE dukklämmare | Backhaus STILLE Towel Clamps

Towel clamps in two different versions. Sharp/sharp or ball/socket for paper towels.

107-10910	Sharp/Sharp	11 cm 4.25"
107-10912	Ball/socket	13.5 cm 5.25"



**Babcock STILLE fattningstång | Babcock STILLE Intestinal Grasping Forceps**

Intestinal grasping forceps in two different lengths.




107-12370

107-12370-16	16 cm 6.25"
107-12370-20	20 cm 8"



**Duval STILLE fattningstång | Duval STILLE Intestinal Grasping Forceps**

Intestinal grasping forceps.



107-12563

107-12563	18 cm 7"
-----------	----------





## Allis STILLE fattningstång | Allis STILLE Intestinal Grasping Forceps

Tissue-grasping forceps for secure grip of tissue.




107-13120

107-13120-15	5x6 teeth	15 cm 6"
107-13120-21	5x6 teeth	21 cm 8.25"



**Allis STILLE fattningstång | Allis STILLE Intestinal Grasping Forceps**

Atraumatic tissue-grasping forceps.




107-13130

107-13130-15 Non-traumatic 15,5 cm 6"


107-13130-19 Non-traumatic 19 cm 7.5"

107-13130-25 Non-traumatic 25 cm 10"



**STILLE-Norell spiraltång | STILLE-Norell IUD Forceps**

IUD forceps with slightly curved, flattened jaws.



107-12491

107-12491 slightly curved 27 cm 10.5"



### Maier STILLE tvättång | Maier STILLE Sponge Holding Forceps

Sponge holding forceps with serrated jaws. Available straight and curved.

107-12478	With ratchet	Straight	26 cm 10.25"
107-12479	With ratchet	Curved	26 cm 10.25"




### Foerster STILLE ringtång | Foerster STILLE Sponge Holding Forceps

Sponge holding forceps with several different loops and lengths.



107-12485



107-12487



107-12486




107-12488

107-12489	Serrated	Straight	18 cm 7"
107-12490	Serrated	Curved	18 cm 7"
107-12485	Serrated	Straight small	24 cm 9.5"
107-12486	Serrated	Curved small	24 cm 9.5"
107-12487	Serrated	Straight large	24 cm 9.5"
107-12488	Serrated	Curved large	24 cm 9.5"



**Schroeder STILLE klotång | Schroeder STILLE  
Tenaculum Forceps**

Uterine tenaculum forceps.



107-13280

107-13280 24 cm 9.5"





**Klotänger** | Tenaculum forceps | Hakenzangen | Pinze per il collo dell'utero | Pinzas para coger el cuello del útero | 子宮支持鉗子

### STILLE uterusmanipulator | STILLE Uterine Manipulator

Designed to give better control of the uterus for improved visualization during laparoscopic procedures. Slightly rounded probe with a drop-shaped tip for minimized risk of tissue trauma.

107-13290 Probe 60 mm Curved 26.5 cm 10.5"

