


Bentänger | Bone cutters and rongeurs | Knochen & hohlmeisselzangen | Cesoie per ossa & pinze sgorbie | Pinzas cortantes & pinzas-gubias | 骨剪刀および骨鉗子(ロンジュール)




Strong joint with screws for easier repair management.

Ergonomic design for less strain in the hand.

Hardness for cutting parts and handles are different from each other.

STILLE plattång | STILLE Flat pliers

General purpose flat pliers with cross-serrated jaws.



105-16015

105-16015 Cross-serrated 17 cm 6.75"



STILLE extraktionstång | STILLE Extraction plier

Original STILLE double-action joint facilitates removal of pins and wires.



105-16020




105-16020-1

105-16020		19 cm 7.5"
105-16020-1	Narrow jaws	19 cm 7.5"




STILLE benavbitare | STILLE Bone Cutting Forceps

Bone cutting forceps with original STILLE double-action joint.



STILLE benavbitare | STILLE Bone Cutting Forceps

Delicate bone cutting forceps with original STILLE double-action joint.



STILLE benavbitare | STILLE Bone Cutting Forceps


Original STILLE double-action angled bone cutting forceps.

105-16811 Curved 19 cm 7.5"



Beyer STILLE gougetång | Beyer STILLE Bone Rongeur

Nasal rongeur with curved, delicate jaws.



105-17030

105-17030 Curved 17 cm 6.75"



Beyer STILLE gougetång | Beyer STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and two different jaw widths.



105-17210



105-17211

105-17210	Curved		18 cm 7"
105-17211	Curved	Delicate	18 cm 7"



Zaufal-Jansen STILLE gougetång | Zaufal-Jansen STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and wide, curved jaws.



105-17215

105-17215	Curved		17 cm 6.75"
-----------	--------	--	-------------



STILLE gougetång | STILLE Bone Rongeur

Bone rongeurs with STILLE double-action joint and two different jaw widths.



105-17216



105-17217

105-17216		Curved	15 cm 6"
105-17217	Delicate jaws	Curved	15 cm 6"



Ruskin-STILLE gougetång | Ruskin-STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and curved jaws.




105-17220

105-17220		Curved	19 cm 7.5"
-----------	--	--------	------------



Luer-STILLE gougetång | Luer-STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and strongly curved wide jaws.




105-17230

105-17230 Curved 22 cm 8,75"



STILLE gougetång | STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and sideways curved jaws.



105-17250

105-17250 Curved 22 cm 8.75"



Leksell STILLE gougetång | Leksell STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and sideways slightly curved jaws.



105-17260




105-17261

105-17260	Curved	Wide jaws	24 cm 9.5"
105-17261	Curved	Narrow jaws	24 cm 9.5"



Leksell STILLE gougetång | Leksell STILLE
Bone Rongeur

Bone rongeur with original STILLE double-action joint and sideways strongly curved jaws.



105-17265


105-17266

105-17265	Curved	Wide jaws	24 cm 9.5"
105-17266	Curved	Narrow jaws	24 cm 9.5"



STILLE gougetång | STILLE Bone Rongeur

Bone rongeur with original STILLE double-action joint and sideways strongly curved jaws in different sizes.



105-17267	Curved jaws	2 mm	24 cm 9.5"
105-17268	Curved jaws	3 mm	24 cm 9.5"
105-17269	Curved jaws	4 mm	24 cm 9.5"



**Semb STILLE benhållartång | Semb STILLE
Bone Holding Forceps**

Bone holding forceps with or without locking mechanism.

105-16150	Without lock	19,5 cm 7.75"
105-16152	With lock	19,5 cm 7.75"



STILLE trådavbitare | STILLE Wire Cutter

Single joint wire cutter with TC-inserts. For hard wire up to 0.4 mm and soft wire up to 0.6 mm.

105-16215 TC inserts 12.5 cm 5"



STILLE trådavbitare | STILLE Wire Cutter

Double joint wire cutter with TC-inserts. For hard wire up to 1.7 mm and soft wire up to 2.0 mm.

105-16223 TC inserts 17.5 cm 7"



STILLE trådavbitare | STILLE Wire Cutter

Strong double joint wire cutter with TC-inserts. For hard wire up to 2.5 mm and soft wire up to 2.8 mm.

105-16225 TC inserts 22 cm 8.75"



Kerrison STILLE laminektomistans | Kerrison STILLE Laminectomy Punch

The popular Kerrison punches with a unique rubber coated handle for improved ergonomics. No spring between the handles improves the accessibility in cleaning and decreases the risk for getting caught with the gloves. All models come with ejector, 130° up cutting and thin footplate.

105-17642-2	18 cm 7"	2 mm
105-17642-3	18 cm 7"	3 mm
105-17642-4	18 cm 7"	4 mm
105-17642-5	18 cm 7"	5 mm
105-17643-2	20 cm 8"	2 mm
105-17643-3	20 cm 8"	3 mm
105-17643-4	20 cm 8"	4 mm
105-17643-5	20 cm 8"	5 mm
105-17644-2	23 cm 9"	2 mm
105-17644-3	23 cm 9"	3 mm
105-17644-4	23 cm 9"	4 mm
105-17644-5	23 cm 9"	5 mm



Kerrison STILLE laminektomistans | Kerrison STILLE Laminectomy Punch

The same fine Kerrison punches are also available with coating for increased durability.

105-17642-2-C	18 cm 7"	2 mm
105-17642-3-C	18 cm 7"	3 mm
105-17642-4-C	18 cm 7"	4 mm
105-17642-5-C	18 cm 7"	5 mm
105-17643-2-C	20 cm 8"	2 mm
105-17643-3-C	20 cm 8"	3 mm
105-17643-4-C	20 cm 8"	4 mm
105-17643-5-C	20 cm 8"	5 mm
105-17644-2-C	23 cm 9"	2 mm
105-17644-3-C	23 cm 9"	3 mm
105-17644-4-C	23 cm 9"	4 mm
105-17644-5-C	23 cm 9"	5 mm



Cushing STILLE bittång | Cushing STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft.
 Available with or without coating. Jaw size 2x10 mm.

105-13921-18	2x10 mm	18 cm 7"	
105-13921-18-C	2x10 mm	18 cm 7"	Coated



Love-Gruenwald STILLE bittång | Love-Gruenwald STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft.
 Available with or without coating. Jaw size 3x10 mm.

105-13931-18	3x10 mm	18 cm 7"	
105-13931-18-C	3x10 mm	18 cm 7"	Coated



Spurling STILLE bittång | Spurling STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Jaw size 4x10 mm.

105-13941-18	4x10 mm	18 cm 7"
105-13941-18-C	4x10 mm	18 cm 7" Coated



Caspar STILLE bittång | Caspar STILLE Rongeur

Rongeur for intervertebral discs, with semi-detachable shaft. Available with or without coating. Fenestrated and serrated jaws in three different sizes.

105-13955-18	2x12 mm	18 cm 7"	
105-13955-18-C	2x12 mm	18 cm 7"	Coated
105-13956-18	3x12 mm	18 cm 7"	
105-13956-18-C	3x12 mm	18 cm 7"	Coated
105-13957-18	4x14 mm	18 cm 7"	
105-13957-18-C	4x14 mm	18 cm 7"	Coated

