



STILLE mikrosax | STILLE Micro Scissors

MIKROSAXAR

Sortimentet består av både raka och böjda modeller i olika längder, bladstorlek och vinklar. Alla modeller är gjorda med STILLE original SuperCut blad för perfekt precision och prestanda i mjuk vävnad. Ledskruven kan tas bort och återmonteras i samband med service. Efter service skall instrumentet fungera med precis samma känsla och pålitlighet som ett nytt instrument.

MICRO SCISSORS

The STILLE Micro scissors consist of both straight and curved models in different lengths, blade sizes and angles. All models are made with the original STILLE Supercut blades for consistent cutting performance in soft tissue. Our joint screw design allows for easy removal during instrument service and repair. After service, the instrument will perform with the same feel and reliability you come to expect from a STILLE instrument.



STILLE mikropincett | STILLE Micro Forceps



STILLE Titan mikropincett | STILLE Titanium Micro Forceps

MIKROPINCETTER


Alla STILLE mikropincetter är gjorda för att ge samma känsla och balans som våra standard-pincetter. Våra duktiga instrumentmakare är fokuserade på de detaljer som ger perfekt fjäderspänning och matchning av spetsarna för att mjukt och enkelt kunna greppa vävnad och suturer. Pincetterna erbjuds i rostfritt stål eller titan, modeller finns med suturplattform, 1x2 klor, ringspetsar och DeBakey atraumatiska kärpincetter. Käftarna kan fås som standard i rostfritt, eller med Diamond Jaws för utökad hållbarhet och greppstyrka.

MICRO FORCEPS

All STILLE Micro Forceps are designed to give the same feel and balance as our standard line of forceps. Our attention to detail by our skilled artisans in Sweden, combines both perfect spring tension and tip alignment for reliable grasping of tissue. They are available in Titanium or stainless steel that include tying platform, 1x2 teeth, ring tips and DeBakey. The tips are available in stainless steel or STILLE Diamond jaws for added durability.



STILLE mikronålförare | STILLE Micro Needle Holder



STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

MIKRONÅLFÖRARE

STILLES mikronålförare är designade för att hålla även de finaste mikrosuturerna. Bara ett minimalt tryck behövs från kirurgens fingrar för ett säkert grepp. Alla modeller är försedda med STILLE Diamond Jaws eller på de kraftigare modellerna hårdmetallplattor. De erbjuds som raka eller böjda, med- eller utan låsning och i ett flertal olika längder. Leddesign Regular för 7-0 suturer och finare eller Streamline för 5-0 och finare (kraftiga modellerna för 2-0 till 4-0).

MICRO NEEDLE HOLDERS

The STILLE Micro Needle Holders are designed to hold the finest sutures in a tight and secure grip. Only minimal pressure will need to be applied by surgeon's fingers. All models are made with STILLE Diamond Jaws and the heavier models with TC-inserts. They can be delivered as straight or curved, with or without ratchet and in several different lengths. Regular box lock for 7-0 sutures and finer, Streamline box lock for 5-0 sutures and finer. (heavy models for 2-0 to 4-0).



STILLE SuperCut mikrosax | STILLE SuperCut Micro Scissors

Fine micro scissors with straight or curved blades with sharp or blunt tips.
All models with **SuperCut** blades.



117-A110
117-A112




117-A111
117-A113

117-A110-15	Sharp tips	Straight	15 cm 6"
117-A110-18	Sharp tips	Straight	18 cm 7"
117-A110-21	Sharp tips	Straight	21 cm 8.25"
117-A112-15	Blunt tips	Straight	15 cm 6 in
117-A112-18	Blunt tips	Straight	18 cm 7 in
117-A112-21	Blunt tips	Straight	21 cm 8.25 in
117-A111-15	Blunt tips	Curved	15 cm 6"
117-A111-18	Blunt tips	Curved	18 cm 7"
117-A111-21	Blunt tips	Curved	21 cm 8.25"
117-A113-15	Sharp tips	Curved	15 cm 6 in
117-A113-18	Sharp tips	Curved	18 cm 7 in
117-A113-21	Sharp tips	Curved	21 cm 8.25 in




STILLE SuperCut mikrokärlsax | STILLE SuperCut Micro Vascular Scissors


SuperCut micro scissors available in five different blade angles
for fine vascular surgery. Blade length 11 mm.




117-A120-25
117-A121-25




117-A120-45
117-A121-45



117-A120-60
117-A121-60



117-A120-90
117-A121-90




117-A120-125
117-A121-125

117-A120-25	Angled 25°	18 cm 7"
117-A120-45	Angled 45°	18 cm 7"
117-A120-60	Angled 60°	18 cm 7"
117-A120-90	Angled 90°	18 cm 7"
117-A120-125	Angled 125°	18 cm 7"
117-A121-25	Angled 25°	21 cm 8.25 in
117-A121-45	Angled 45°	21 cm 8.25 in
117-A121-60	Angled 60°	21 cm 8.25 in
117-A121-90	Angled 90°	21 cm 8.25 in
117-A121-125	Angled 125°	21 cm 8.25 in

STILLE SuperCut mikrokärlsax | STILLE SuperCut Micro Vascular Scissors

Extra fine blade version of STILLE **SuperCut** vascular micro scissors.
Blade length 7.5 mm.



117-A125-25	Angled 25°	18 cm 7"
117-A125-45	Angled 45°	18 cm 7"
117-A125-60	Angled 60°	18 cm 7"
117-A125-90	Angled 90°	18 cm 7"
117-A125-125	Angled 125°	18 cm 7"



STILLE mikropincett | STILLE Micro Forceps

Well balanced and ergonomic micro forceps of ultra high quality. With tying platform and straight or curved jaws. 0.4 mm tips.




117-B110	117-B111		
117-B110-15	Straight	15 cm 6"	
117-B110-18	Straight	18 cm 7"	
117-B110-21	Straight	21 cm 8.25"	
117-B110-23	Straight	23 cm 9 in	
117-B111-15	Curved	15 cm 6"	
117-B111-18	Curved	18 cm 7"	
117-B111-21	Curved	21 cm 8.25"	



STILLE mikropincett | STILLE Micro Forceps

Same fine micro forceps, with Diamond Jaws, for increased tip durability and secure grip of fine needles. 0.4 mm tips.



117-B120	117-B121		
117-B120-15	Diamond Jaws	Straight	15 cm 6"
117-B120-18	Diamond Jaws	Straight	18 cm 7"
117-B120-21	Diamond Jaws	Straight	21 cm 8.25"
117-B120-23	Diamond Jaws	Straight	23 cm 9 in
117-B121-15	Diamond Jaws	Curved	15 cm 6"
117-B121-18	Diamond Jaws	Curved	18 cm 7"
117-B121-21	Diamond Jaws	Curved	21 cm 8.25"



STILLE mikropincett | STILLE Micro Forceps

Micro forceps with Diamond Jaws and Counter-balance. 0.4 mm tips.



117-B130	117-B131		
117-B130-18	Diamond Jaws	Straight	18 cm 7"
117-B130-21	Diamond Jaws	Straight	21 cm 8.25"
117-B130-23	Diamond Jaws	Straight	23 cm 9 in
117-B131-18	Diamond Jaws	Curved	18 cm 7"
117-B131-21	Diamond Jaws	Curved	21 cm 8.25"



STILLE mikro ringpincett | STILLE Micro Ring Forceps

Counter-balanced micro forceps with fine Diamond Jaws ring tips.



117-B140

117-B140-18	1.0x0.5 mm	Diamond Jaws	Straight	18 cm 7"
117-B140-21	1.0x0.5 mm	Diamond Jaws	Straight	21 cm 8.25"
117-B140-23	1.0x0.5 mm	Diamond Jaws	Straight	23 cm 9 in



STILLE mikro klopincett | STILLE Micro Tissue Forceps

Round handle, with 1x2 fine teeth on 0.4 mm tips.

117-B150-15	1x2 fine teeth	Straight	15 cm 6"
117-B150-18	1x2 fine teeth	Straight	18 cm 7"
117-B150-21	1x2 fine teeth	Straight	21 cm 8.25"



STILLE DeBakey mikropincett | STILLE Micro DeBakey Forceps

Round handle, with 1.2 mm fine atraumatic DeBakey tips.

117-B160-18	1.2 mm	DeBakey tips	Straight	18 cm 7 in
117-B160-21	1.2 mm	DeBakey tips	Straight	21 cm 8.25 in



STILLE DeBakey mikropincett | STILLE Micro DeBakey Forceps

Round handle, with 1.2 mm fine atraumatic DeBakey tips, counter-balanced.

117-B165-18	1.2 mm	DeBakey tips	Straight	18 cm 7 in
117-B165-21	1.2 mm	DeBakey tips	Straight	21 cm 8.25 in



STILLE Titan mikropincett | STILLE Titanium Micro Forceps

Round handle, with tying platform, straight or curved 0.4 mm tips.

117-B120-18-TI	Diamond Jaws	Straight	18 cm 7 in
117-B120-21-TI	Diamond Jaws	Straight	21 cm 8.25 in
117-B120-23-TI	Diamond Jaws	Straight	23 cm 9 in
117-B120-25-TI	Diamond Jaws	Straight	25 cm 10 in
117-B121-18-TI	Diamond Jaws	Curved	18 cm 7 in
117-B121-21-TI	Diamond Jaws	Curved	21 cm 8.25 in
117-B121-23-TI	Diamond Jaws	Curved	23 cm 9 in
117-B121-25-TI	Diamond Jaws	Curved	25 cm 10 in



STILLE Titan mikropincett | STILLE Titanium Micro Forceps

Round handle, with tying platform, counter-balanced, straight or curved 0.4 mm tips.

117-B130-18-TI	Diamond Jaws	Straight	18 cm 7 in
117-B130-21-TI	Diamond Jaws	Straight	21 cm 8.25 in
117-B131-18-TI	Diamond Jaws	Curved	18 cm 7 in
117-B131-21-TI	Diamond Jaws	Curved	21 cm 8.25 in



NEW
Titanium

Mikropincetter | Micro forceps | Mikropinzetten | Micro pinze | Micro pinzas | マイクロ鑷子



STILLE Titan mikro ringpincett | STILLE Titanium Micro Ring Forceps

Round handle, with fine ring tips, counter-balanced.

117-B140-18-TI	Diamond Jaws	1.0x0.5 mm ring tips	Straight	18 cm 7 in
117-B140-21-TI	Diamond Jaws	1.0x0.5 mm ring tips	Straight	21 cm 8.25 in
117-B140-23-TI	Diamond Jaws	1.0x0.5 mm ring tips	Straight	23 cm 9 in



STILLE Titan DeBakey mikropincett | STILLE Titanium Micro DeBakey Forceps

Round handle, with 1.2 mm fine atraumatic DeBakey tips.

117-B160-18-TI	1.2 mm DeBakey tips	Straight	18 cm 7 in
117-B160-21-TI	1.2 mm DeBakey tips	Straight	21 cm 8.25 in





STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock** for sutures 7-0 and smaller. Without ratchet.



117-C110-15	Diamond Jaws	Straight	15 cm 6"
117-C110-18	Diamond Jaws	Straight	18 cm 7"
117-C110-21	Diamond Jaws	Straight	21 cm 8.25"
117-C110-23	Diamond Jaws	Straight	23 cm 9"
117-C111-15	Diamond Jaws	Curved	15 cm 6"
117-C111-18	Diamond Jaws	Curved	18 cm 7"
117-C111-21	Diamond Jaws	Curved	21 cm 8.25"
117-C111-23	Diamond Jaws	Curved	23 cm 9"

STILLE mikro nålförare | STILLE Micro Needle Holder


Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock** for sutures 5-0 and smaller. Without ratchet.

117-C115-15	Diamond Jaws	Without Ratchet	Straight	15 cm 6 in
117-C115-18	Diamond Jaws	Without Ratchet	Straight	18 cm 7 in
117-C115-21	Diamond Jaws	Without Ratchet	Straight	21 cm 8.25 in
117-C115-23	Diamond Jaws	Without Ratchet	Straight	23 cm 9 in
117-C116-15	Diamond Jaws	Without Ratchet	Curved	15 cm 6 in
117-C116-18	Diamond Jaws	Without Ratchet	Curved	18 cm 7 in
117-C116-21	Diamond Jaws	Without Ratchet	Curved	21 cm 8.25 in
117-C116-23	Diamond Jaws	Without Ratchet	Curved	23 cm 9 in



STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock** for sutures 7-0 and smaller. With ratchet.



117-C120-15	Diamond Jaws	Straight	15 cm 6"
117-C120-18	Diamond Jaws	Straight	18 cm 7"
117-C120-21	Diamond Jaws	Straight	21 cm 8.25"
117-C120-23	Diamond Jaws	Straight	23 cm 9"
117-C121-15	Diamond Jaws	Curved	15 cm 6"
117-C121-18	Diamond Jaws	Curved	18 cm 7"
117-C121-21	Diamond Jaws	Curved	21 cm 8.25"
117-C121-23	Diamond Jaws	Curved	23 cm 9"



STILLE mikro nålförare | STILLE Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock** for sutures 5-0 and smaller. With ratchet.

117-C125-15	Diamond Jaws	With Ratchet	Straight	15 cm 6 in
117-C125-18	Diamond Jaws	With Ratchet	Straight	18 cm 7 in
117-C125-21	Diamond Jaws	With Ratchet	Straight	21 cm 8.25 in
117-C125-23	Diamond Jaws	With Ratchet	Straight	23 cm 9 in
117-C126-15	Diamond Jaws	With Ratchet	Curved	15 cm 6 in
117-C126-18	Diamond Jaws	With Ratchet	Curved	18 cm 7 in
117-C126-21	Diamond Jaws	With Ratchet	Curved	21 cm 8.25 in
117-C126-23	Diamond Jaws	With Ratchet	Curved	23 cm 9 in



STILLE mikro nålförare | STILLE Micro Needle Holder

Extra heavy model for handling of larger needles with TC-inserts. Intended for 2/0 suture and finer.



117-C130
117-C140

117-C130-19	TC-inserts	Without lock	19 cm 7.5"
117-C130-23	TC-inserts	Without lock	23 cm 9"
117-C130-25	TC-inserts	Without lock	25 cm 9.75"
117-C140-19	TC-inserts	With lock	19 cm 7.5"
117-C140-23	TC-inserts	With lock	23 cm 9"
117-C140-25	TC-inserts	With lock	25 cm 9.75"



STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

Round handle, fine 0.8 mm Diamond Jaws and **Regular Box Lock**
for sutures 7-0 and smaller. With or without ratchet.

117-C110-18-TI	Diamond Jaws	Without Ratchet	Straight	18 cm 7 in
117-C110-21-TI	Diamond Jaws	Without Ratchet	Straight	21 cm 8.25 in
117-C110-23-TI	Diamond Jaws	Without Ratchet	Straight	23 cm 9 in
117-C111-15-TI	Diamond Jaws	Without Ratchet	Curved	15 cm 6 in
117-C111-18-TI	Diamond Jaws	Without Ratchet	Curved	18 cm 7 in
117-C111-21-TI	Diamond Jaws	Without Ratchet	Curved	21 cm 8.25 in
117-C120-18-TI	Diamond Jaws	With Ratchet	Straight	18 cm 7 in
117-C120-21-TI	Diamond Jaws	With Ratchet	Straight	21 cm 8.25 in
117-C120-23-TI	Diamond Jaws	With Ratchet	Straight	23 cm 9 in
117-C121-18-TI	Diamond Jaws	With Ratchet	Curved	18 cm 7 in
117-C121-21-TI	Diamond Jaws	With Ratchet	Curved	21 cm 8.25 in
117-C121-23-TI	Diamond Jaws	With Ratchet	Curved	23 cm 9 in



NEW
Titanium

Mikronålförare | Micro needle holders | Mikronadelhalter | Micro portaghi |
Micro porta-agujas | マイクロ持針器



STILLE Titan mikronålförare | STILLE Titanium Micro Needle Holder

Round handle, 1.2 mm Diamond Jaws and **Streamline Box Lock**
for sutures 5-0 and smaller. With or without ratchet.

117-C115-18-TI	Diamond Jaws	Without Ratchet	Straight	18 cm 7 in
117-C115-21-TI	Diamond Jaws	Without Ratchet	Straight	21 cm 8.25 in
117-C115-23-TI	Diamond Jaws	Without Ratchet	Straight	23 cm 9 in
117-C125-18-TI	Diamond Jaws	With Ratchet	Straight	18 cm 7 in
117-C125-21-TI	Diamond Jaws	With Ratchet	Straight	21 cm 8.25 in
117-C125-23-TI	Diamond Jaws	With Ratchet	Straight	23 cm 9 in
117-C126-18-TI	Diamond Jaws	With Ratchet	Curved	18 cm 7 in
117-C126-21-TI	Diamond Jaws	With Ratchet	Curved	21 cm 8.25 in
117-C126-23-TI	Diamond Jaws	With Ratchet	Curved	23 cm 9 in

