


Pincetter | Forceps | Pinzetten | Pinze | Pinzas |


鑷子 (ピンセット)





Adson STILLE pincett | Adson STILLE Tissue Forceps

Tissue forceps with very delicate points and 1x2 fine teeth. Used in most specialties for delicate work.



102-9019


102-9022
102-9022-15

102-9019	Very delicate	1x2 teeth	12 cm 4.75"
102-9022	Cross-serrated	1x2 teeth	12 cm 4.75"
102-9022-15	Cross-serrated	1x2 teeth	15 cm 6"



Adson-Callison STILLE pincett | Adson-Callison STILLE Tissue Forceps

Tissue forceps with tying platform and 1x2 fine teeth.



102-9021

102-9021-1

102-9021	Tying platform	1x2 teeth	12 cm 4.75"
102-9021-1	Delicate tying platform	1x2 teeth	12 cm 4.75"



Adson-Brown STILLE pincett | Adson-Brown STILLE Tissue Forceps

Tissue forceps with side-grasping teeth for secure, atraumatic grasping of soft tissue or cartilage.



102-9023-12

102-9023-12	Side-grasping teeth	12 cm 4.75"
-------------	---------------------	-------------

Gillies STILLE pincett | Gillies STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth and stabilizing pin. An excellent forceps for the surgeon who wants an Adson-size tip, but with greater length.



102-9027

102-9027	Serrated	1x2 fine teeth	15.5 cm 6.25"
----------	----------	----------------	---------------



STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth and stabilizing pin. Delicate model.



102-9028

102-9028	Serrated	1x2 fine teeth	15.5 cm 6.25"
----------	----------	----------------	---------------



STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth.



102-9050-13

102-9050-13	1x2 fine teeth	13 cm 5"
-------------	----------------	----------

STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 1x2 teeth.



102-9060

102-9060-15	1x2 teeth	15 cm 6"
102-9060-17	1x2 teeth	17 cm 6.75"
102-9060-20	1x2 teeth	20 cm 8"



STILLE pincett | STILLE Tissue Forceps

Tissue forceps with 2x3 teeth. Used when more grasping power than 1x2 teeth is needed.



102-9062

102-9062-15	2x3 teeth	15 cm 6"
102-9062-17	2x3 teeth	17 cm 6.75"
102-9062-20	2x3 teeth	20 cm 8"





Adlercreutz STILLE pincett | Adlercreutz STILLE Tissue Forceps

Tissue forceps with 4x5 fine teeth. Very useful in neurosurgery to hold the dura, and in cardiovascular surgery to hold the pericardium.



102-9248-15	4x5 fine teeth	15 cm 6"
102-9248-20	4x5 fine teeth	20 cm 8"
102-9248-25	4x5 fine teeth	25 cm 10"

Waugh-Brophy STILLE pincett | Waugh-Brophy STILLE Tissue Forceps

Tissue forceps with 1x2 fine teeth. Very fine tips, like Adson. Useful when long and slender instruments are needed.



102-9252

102-9252-20	1x2 teeth	Cross-serrated	20 cm 8"
102-9252-25	1x2 teeth	Cross-serrated	25 cm 10"
102-9252-30	1x2 teeth	Cross-serrated	30 cm 12"



Russian STILLE pincett | Russian STILLE Tissue Forceps

Forceps for secure grip of tissue. Used when a better grip than with a dressing forceps is needed, and when damage to the tissue must be avoided. For hard and slippery tissue.



102-9270-20


102-9270-20	20 cm 8"
-------------	----------





Adson STILLE anatomisk pincett | Adson STILLE Dressing Forceps


General-purpose dressing forceps with various types of jaw surfaces.



102-9020	Serrated	12 cm 4.75"
102-9020-15	Serrated	15 cm 6"
102-9020-1	Diamond jaws	12 cm 4.75"
102-9024	TC jaws cross-serrated	12 cm 4.75"

STILLE anatomisk pincett | STILLE Dressing Forceps

Dressing forceps with square shaped and serrated points used in all surgical specialties.




102-9030-13	Square points	13 cm 5"
102-9030-15	Square points	15 cm 6"
102-9030-17	Square points	17 cm 6.75"
102-9030-20	Square points	20 cm 8"
102-9030-25	Square points	25 cm 10"



McIndoe STILLE anatomisk pincett | McIndoe STILLE Dressing Forceps

Dressing forceps with stabilizing pin and delicate jaws for grasping fine tissue.



102-9025	Serrated	15.5 cm 6.25"
102-9025-1	Diamond Jaws	15.5 cm 6.25"

STILLE anatomisk pincett | STILLE Dressing Forceps

Dressing forceps with rounded points and serrated jaws.





102-9040-15	Rounded points	15 cm 6"
102-9040-17	Rounded points	17 cm 6.75"
102-9040-20	Rounded points	20 cm 8"
102-9040-25	Rounded points	25 cm 10"
102-9040-30	Rounded points	30 cm 12"





**Crafoord-STILLE anatomisk pincett |
Crafoord-STILLE Dressing Forceps**

Dressing forceps with rounded points. The perfect balance and delicate grip gives this instrument a unique feeling.





102-9077-20	Serrated	20 cm 8"
102-9078-20	Diamond jaws	20 cm 8"



**Waugh-Brophy STILLE anatomisk pincett |
Waugh-Brophy STILLE Dressing Forceps**

Dressing forceps with very fine points and two different jaw surfaces. The cross-serration gives excellent holding power.



102-9254	Cross- serrated	20 cm 8"
102-9255-20	Diamond jaws	20 cm 8"
102-9255-25	Diamond jaws	25 cm 10"

Olivecrona-STILLE durapincett | Olivecrona-STILLE Dura Forceps

Dura forceps with 2x3 fine teeth. The nonpenetrating tips give the surgeon maximum grip without damaging the tissue.




102-9244-15	2x3 fine teeth	15 cm 6"
102-9244-20	2x3 fine teeth	20 cm 8"
102-9244-25	2x3 fine teeth	25 cm 10"



STILLE irispincett | STILLE Iris Forceps

Short delicate forceps with either 1x2 fine teeth or serrated jaws.



102-9460




102-9470

102-9460	Serrated	10 cm 4"
102-9470	1x2 teeth	10 cm 4"



STILLE De Bakey kärlpincett | STILLE De Bakey Vascular Forceps

Atraumatic forceps with De Bakey jaws for handling of blood vessels and other delicate tissue. Can be insulated on request.




102-9300 102-9302 102-9312 102-9314 102-9320
102-9304 102-9316 102-9318 102-9322

102-9340	De Bakey type	1 mm	15 cm 6"
102-9342	De Bakey type	1 mm	20 cm 8"
102-9354	De Bakey type	1.5 mm	25 cm 10"
102-9300	De Bakey type	2 mm	15 cm 6"
102-9302	De Bakey type	2 mm	20 cm 8"
102-9304	De Bakey type	2 mm	25 cm 10"
102-9312	De Bakey type	2.7 mm	15 cm 6"
102-9314	De Bakey type	2.7 mm	20 cm 8"
102-9316	De Bakey type	2.7 mm	25 cm 10"
102-9318	De Bakey type	2.7 mm	30 cm 11.75"
102-9320	De Bakey type	3.2 mm	20 cm 8"
102-9322	De Bakey type	3.2 mm	25 cm 10"



STILLE De Bakey kärlpincett | STILLE De Bakey Vascular Forceps

Angled atraumatic forceps with De Bakey jaws for delicate tissue.




102-9308 De Bakey type 2 mm 20 cm 8"



STILLE De Bakey kärlpincett | STILLE De Bakey Vascular Forceps

Atraumatic forceps with De Bakey jaws for delicate tissue, non serrated handles. Stop pin for preventing from gripping the tissue too hard.



102-9330

102-9330 De Bakey type 2 mm 20 cm 8"





**Waugh-Brophy STILLE isolerad diatermipincett |
 Waugh-Brophy STILLE insulated Diathermy
 Forceps**

Insulated version of the popular Waugh-Brophy forceps. Used in monopolar diathermy with 1x2 teeth.



102-9852

102-9852-20	1x2 teeth	Cross-serrated	20 cm 8"
102-9852-25	1x2 teeth	Cross-serrated	25 cm 10"
102-9852-30	1x2 teeth	Cross-serrated	30 cm 12"

**Björk-STILLE diatermipincetter | Björk-STILLE
 Diathermy Forceps**

Insulated forceps for monopolar diathermy. Serrations allow grasping of delicate structures.



102-9845

102-9845-15	Insulated	15 cm 6"
102-9845-20	Insulated	20 cm 8"
102-9845-25	Insulated	25cm 10"

**Waugh-Brophy STILLE isolerad diatermipincett |
 Waugh-Brophy STILLE insulated Diathermy
 Forceps**

Insulated version of the popular Waugh-Brophy forceps. Used in monopolar diathermy. With Diamond Jaws.



102-9855

102-9855-20	Diamond jaws	20 cm 8"
102-9855-25	Diamond jaws	25 cm 10"



Björk–STILLE knutpincett | Björk–STILLE Knot-Holding Forceps

Same as Björk–STILLE diathermy forceps, but without insulation.



102-9846-20

