


Saxar | Scissors | Scheren | Forbici | Tijeras | 剪刀 (はさみ)



Smooth finger rings for comfort.

Shanks designed for feeling and balance.

Blades tempered to the correct hardness.

Strong joint for correct tension between blades.

Extreme precision right down to the tip of the edges.



STILLE Ögonsax | STILLE Eye Scissors

ORIGINAL

Saxar har alltid tillhört STILLES storsäljare alltsedan de började tillverkas på 1800-talet. Genom åren har de ständigt förbättrats och fler specialmodeller har tillkommit. Idag finns de i totalt ca 200 varianter i alla storlekar och utformningar.

Originalsaxarna tillverkas i rostfritt stål av högsta kvalitet med högglosspolerad yta. De är välbalanserade med god skärstyrka ända ut i spetsarna och ligger bra i handen. Det finns originalsaxar för alla typer av mjukvävnad och de kan även användas för att klippa gasbinda, tråd och suturer.

SUPERCUT

Saxarnas sax uppfanns 1982. Många kirurger anser att SuperCut är världens absolut bästa kirurgisax. Den passar alla typer av mjuk vävnad och precisionen och prestandan är välkänd bland kirurger över hela världen.

Det revolutionerande med SuperCut är skärytorna som samarbetar på ett unikt sätt. Det ena bladet är saxslipat och det andra är knivslipat så att saxen skär sig fram genom vävnaden. Det ger maximal skärstyrka med minimal kraft och ett perfekt snitt. Slipningen gör saxen lämplig för både höger- och vänsterhänta. Trots många år på marknaden har ingen lyckats skapa en kopia som kan överträffa originalet.

STILLES nya beläggning på SuperCut görs med ett mycket tunt diamantliknande skikt av s.k. DLC-coating, som reducerar risken för missfärgning, ökar instrumentets hållbarhet och bidrar till en bättre miljö.

ORIGINAL

Scissors have always been STILLE's bestseller ever since they were first produced in the 1800's. Through the years they have been continually improved and more special models have been introduced. Currently, there are in total about 200 different varieties in all sizes and shapes.

Our original scissors are manufactured of the highest quality, polished stainless steel. They are well balanced with outstanding cutting power right down to the tips and they fit well in the surgeon's hand. There are original scissors for all types of soft tissue and they can also be used to cut gauze, wire and sutures.

SUPERCUT

These unique scissors were introduced in 1982. Many surgeons believe that the SuperCut are unsurpassed as the world's best surgical scissors. They are suitable for all types of soft tissue and the precision and performance are well known among surgeons worldwide.

The revolutionary SuperCut scissors feature cutting surfaces that work together in a unique way. One blade is a scissor blade and the other is a knife blade to enable the scissors to cut precisely through the soft tissue. They provide maximum cutting power with minimal force and produce a perfect cut. The sharpening makes the scissors suitable for both right and left handed users. Despite many years on the market, other brands have not been able to create an imitation that can surpass the original SuperCut scissors.

The new STILLE coating process is based on a thin layer diamond-like-carbon proprietary technology, aiming to reduce discoloration, improve instrument durability and help to create a more sustainable environment.



STILLE-Mayo SuperCut XE dissektionssax |
STILLE-Mayo SuperCut XE Dissecting Scissors

DIAMOND SUPERCUT

Våra nya STILLE Diamond SuperCut saxar tar kirurgiska saxar till en ny nivå. Vi välkomnar er att uppleva dissektion, klippegenskaper och hållbarhet som ni aldrig tidigare upplevt. Diamond SuperCut är handgjorda av våra skickliga instrumentmakare i Eskilstuna. Saxarna är helbelagda med en ny diamentliknande beläggning, som ökar hållbarheten, minskar behovet av återkommande slipning, är mer motståndskraftig mot korrosion och är mer miljövänlig. Den nya Diamond-beläggningen är validerad för att tåla de allt vanligare rengöringsprocesserna med väteperoxid.



Metzenbaum STILLE Diamond SuperCut dissektionssax
| Metzenbaum STILLE Diamond SuperCut Dissecting Scissors

ERGONOMISKA SAXAR

En naturlig fortsättning på arbetet med att skapa perfekta saxar har varit att vidareutveckla deras ergonomiska egenskaper. Kirurgiska saxar hålls normalt i ett precisionsgrepp med tummen och lång- eller ringfingret. Vid mer tidskrävande ingrepp tröttnas handen ut snabbare.

Lösningen är att vinkla skänklarna så att handen kan arbeta närmare sitt normala viloläge. Med ett större mellanrum mellan fingeröglorna motverkas spänningar och trötthet och handen orkar arbeta längre. Samtliga ergonomiska saxar är av modellen SuperCut med sina unika egenskaper. Sax- och knivslipade blad samverkar för maximal skärstyrka med minimal kraft och exakt precision.


DIAMOND SUPERCUT

The new STILLE Diamond SuperCut redefines the industry standards and breaks the limits in instrument performance. We welcome you to experience dissection, cutting performance and durability like you've never felt before. The Diamond SuperCut series of scissors are hand-made to perfection in Sweden. The scissors are fully coated with a new proprietary diamond-like-layer. It not only increases the life expectancy, it reduces the need of sharpening, is more resistant to corrosion and it's more environmental friendly. The new Diamond coating is validated to withstand the more common hydro peroxide processes.

ERGONOMIC SCISSORS

A natural extension of the ongoing effort to create the perfect scissors has been the further development of their ergonomic properties. Surgical scissors are normally held in a precise grip using the thumb and middle or ring finger. For more time-consuming procedures this can cause the hand to become fatigued faster.


The solution is to angle the shanks so that the hand can work closer to its normal resting position. With a larger gap between the finger loops, tension and fatigue is countered and the hand can work more comfortably for a longer period. All ergonomic scissors are SuperCut models with their unique characteristics. Scissor and knife sharpened blades work together for maximum cutting power with minimal effort and exact precision.



Stevens-STILLE SuperCut XE ergonomisk tenotomisax | Stevens-STILLE SuperCut XE Ergonomic Tenotomy Scissors

STILLE ögonsax | STILLE Eye Scissors



Delicate eye scissors with sharp or blunt tips and large finger rings. These are also called iris scissors and are used for nerves, tendons and vessels. They are intended for facial surgery and any type of delicate surgery. Beveled blades.



101-8001	Straight	Sharp	9.5 cm 3.75"
101-8002	Straight	Blunt	9.5 cm 3.75"
101-8003	Curved	Sharp	9.5 cm 3.75"
101-8004	Curved	Blunt	9.5 cm 3.75"

STILLE ögonsax | STILLE Eye Scissors


Eye scissors for delicate work. Large finger rings for a comfortable grip. Beveled blades.



101-8010	Straight	Sharp	10 cm 4"
101-8011	Straight	Blunt	10 cm 4"
101-8015	Curved	Sharp	10 cm 4"
101-8016	Curved	Blunt	10 cm 4"

STILLE SuperCut XE ögonsax | STILLE SuperCut XE Eye Scissors

SuperCut XE eye scissors with sharp or blunt tips and extremely sharp edges. Beveled blades.




101-8010-1	Straight	Sharp	10 cm 4"	
101-8011-1	Straight	Blunt	10 cm 4"	
101-8015-1	Curved	Sharp	10 cm 4"	
101-8015-1-C	Curved	Sharp	10 cm 4"	Diamond
101-8016-1	Curved	Blunt	10 cm 4"	



STILLE SuperCut XE ögonlockssax | STILLE SuperCut XE Blepharoplasty Scissors

Short delicate SuperCut XE scissors suitable for blepharoplasty.



101-8018-1	Straight	9.5 cm 3.75"
101-8019-1	Curved	9.5 cm 3.75"




STILLE SuperCut XE ergonomisk ögonlockssax | STILLE SuperCut XE Ergonomic Blepharoplasty Scissors

Blepharoplasty SuperCut XE scissors with ergonomically designed shanks for a more relaxed grip.



101-8019-2	Curved	9.5 cm 3.75"
------------	--------	--------------



Ragnell STILLE SuperCut XE ögonlockssax | Ragnell STILLE SuperCut XE Blepharoplasty Scissors

Short delicate SuperCut XE scissors with blunt, flat tips for blepharoplasty.



101-8031-1	Curved	9.5 cm 3.75"
------------	--------	--------------




Ragnell STILLE SuperCut XE ergonomisk ögonlockssax | Ragnell STILLE SuperCut XE Ergonomic Blepharoplasty Scissors

Blepharoplasty SuperCut XE scissors with ergonomically designed shanks for a more relaxed grip.



101-8031-2	Curved	9.5 cm 3.75"
------------	--------	--------------



Kaye STILLE dissektionsax | Kaye STILLE Dissecting Scissors

Scissors with serrated blades for "holding" the tissue, making it easier to cut exactly as far as intended.



101-8020




101-8022



101-8021

101-8020	Blunt serrated	Curved	10 cm 4"
101-8021	Blunt serrated	Curved	12.5 cm 5"
101-8022	Sharp serrated	Curved	12.5 cm 5"



Kaye STILLE SuperCut XE dissektionsax | Kaye STILLE SuperCut XE Dissecting Scissors

SuperCut XE version for improved cutting performance.



101-8021-1

101-8021-1	Blunt serrated	Curved	12.5 cm 5"
------------	----------------	--------	------------




Jabaley STILLE dissektionsax | Jabaley STILLE Dissecting Scissors

Sharp scissors suitable for hand surgery, where very strong scissors are used. Beveled blades.




101-8060



101-8062


101-8060	Straight	12.5 cm 5"
101-8062	Curved	12.5 cm 5"





Jabaley STILLE SuperCut XE dissektionssax | Jabaley STILLE SuperCut XE Dissecting Scissors


SuperCut XE scissors for dissection of tough tissue in hand surgery. Beveled blades.



101-8060-1


101-8062-1

101-8060-1	Straight	12.5 cm 5"	
101-8060-1-C	Straight	12.5 cm 5"	Diamond
101-8062-1	Curved	12.5 cm 5"	



Jabaley STILLE SuperCut XE ergonomisk dissektionssax | Jabaley STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergonomic model to avoid fatigue in the hand. Beveled blades.



101-8060-2


101-8062-2

101-8060-2	Straight	12.5 cm 5"	
101-8062-2	Curved	12.5 cm 5"	



Mayo STILLE dissektionssax | Mayo STILLE Dissecting Scissors

The most used scissors, for all types of surgery. Used for routine dissection and cutting in tissue. Beveled blades.



101-8080-15

101-8082-15


101-8080-17

101-8082-17

101-8080-15	Straight	15 cm 6"	
101-8080-17	Straight	17 cm 6.75"	
101-8082-15	Curved	15 cm 6"	
101-8082-17	Curved	17 cm 6.75"	

Mayo STILLE SuperCut XE dissektionssax | Mayo STILLE SuperCut XE Dissecting Scissors



SuperCut XE dissecting scissors for tough tissue. Beveled blades.



101-8081-15	Straight	15 cm 6"	
101-8083-15	Curved	15 cm 6"	
101-8083-15-C	Curved	15 cm 6"	Diamond

STILLE-Mayo dissektionssax | STILLE-Mayo Dissecting Scissors


Sturdy dissecting scissors with rounded blades.




101-8084-15	Rounded blades	Straight	15 cm 6"
101-8084-17	Rounded blades	Straight	17 cm 6.75"
101-8084-19	Rounded blades	Straight	19 cm 7.5"
101-8086-15	Rounded blades	Curved	15 cm 6"
101-8086-17	Rounded blades	Curved	17 cm 6.75"
101-8086-19	Rounded blades	Curved	19 cm 7.5"

STILLE-Mayo SuperCut XE dissektionssax | STILLE-Mayo SuperCut XE Dissecting Scissors

Sturdy Super Cut scissors with rounded blades for dissection of tough tissue.




101-8085-15	Rounded blades	Straight	15 cm 6"	
101-8085-17	Rounded blades	Straight	17 cm 6.75"	
101-8085-17-C	Rounded blades	Straight	17 cm 6.75"	Diamond
101-8085-19	Rounded blades	Straight	19 cm 7.5"	
101-8087-15	Rounded blades	Curved	15 cm 6"	
101-8087-17	Rounded blades	Curved	17 cm 6.75"	
101-8087-17-C	Rounded blades	Curved	17 cm 6.75"	Diamond
101-8087-19	Rounded blades	Curved	19 cm 7.5"	



Mayo-Noble STILLE dissektionssax | Mayo-Noble STILLE Dissecting Scissors

Short and strong dissecting scissors mainly for orthopaedic surgeons. Beveled blades.




101-8090	Straight	18 cm 7"
101-8092	Curved	18 cm 7"



Baby Mayo STILLE dissektionssax | Baby Mayo STILLE Dissecting Scissors

Delicate version of the regular Mayo scissors. Mainly for plastic- and hand surgery. Beveled blades.



101-8100-15	Straight	15 cm 6"
101-8102-15	Curved	15 cm 6"






**Sims STILLE dissektionssax |
Sims STILLE Dissecting Scissors**

Mainly used in gynecological and other pelvic surgery. Beveled blades.



101-8138

101-8138	Curved	20 cm 8"
----------	--------	----------



**Baby Metzenbaum STILLE dissektionssax |
Baby Metzenbaum STILLE Dissecting Scissors**

Short, delicate dissecting scissors with rounded blades. Used in delicate surgery requiring small instruments.



101-8160




101-8162

101-8160	Straight	12,5 cm 5"
101-8162	Curved	12,5 cm 5"

Baby Metzenbaum STILLE SuperCut XE dissektionssax | Baby Metzenbaum STILLE SuperCut XE Dissecting Scissors


Super Cut dissecting scissors for tough tissue. Used in delicate surgery requiring small instruments.



101-8161

101-8163

101-8161	Straight	12,5 cm 5"
101-8163	Curved	12,5 cm 5"



Baby Metzenbaum STILLE Diamond SuperCut dissektionssax | Baby Metzenbaum STILLE Diamond SuperCut Dissecting Scissors

With Diamond coating for increased durability.



101-8163-C

101-8163-C	Curved	12,5 cm 5"
------------	--------	------------



Baby Metzenbaum STILLE SuperCut XE ergonomisk dissektionssax | Baby Metzenbaum STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergonomic model for better control. Rounded blades.




101-8163-2

101-8163-2	Curved	12,5 cm 5"
------------	--------	------------



Metzenbaum STILLE gracil dissektions sax | Metzenbaum STILLE Delicate Dissecting Scissors

Short, delicate dissecting scissors with rounded blades. Used in delicate surgery requiring small instruments.



101-8172

101-8172-18 Curved delicate 18 cm 7"

101-8172-23 Curved delicate 23 cm 9"



Metzenbaum STILLE gracil SuperCut XE dissektions sax | Metzenbaum STILLE Delicate SuperCut XE Dissecting Scissors

Very delicate version with Super Cut edges and fine rounded tips.



101-8173-15



101-8173-18
101-8173-20

101-8173-15 Curved delicate 15 cm 6"

101-8173-18 Curved delicate 18 cm 7"

101-8173-20 Curved delicate 20 cm 8"



Metzenbaum STILLE gracil Diamond SuperCut dissektionssax | Metzenbaum STILLE Delicate Diamond SuperCut Dissecting Scissors

With Diamond coating for increased durability.



101-8173-15-C


101-8173-18-C

101-8173-15-C	Curved delicate	15 cm 6"
101-8173-18-C	Curved delicate	18 cm 7"



Metzenbaum STILLE dissektionssax | Metzenbaum STILLE Dissecting Scissors

Dissecting scissors with rounded blades and tips. Used in all specialties for fine dissection. Metzenbaum is one of STILLE's most popular scissors model and is available in a variety of shapes and sizes.



101-8174-15, 101-8174-18
101-8174-20, 101-8174-23


101-8174-26

101-8174-15	Curved	15 cm 6"
101-8174-18	Curved	18 cm 7"
101-8174-20	Curved	20 cm 8"
101-8174-23	Curved	23 cm 9"
101-8174-26	Curved	26 cm 10.25"



STILLE-Ekberg Metzenbaum gracil SuperCut ergonomisk dissektionssax | STILLE-Ekberg Metzenbaum Delicate SuperCut Ergonomic Dissecting Scissors

Ergonomic STILLE-Ekberg version. Available with blades curved up or down, suitable for both left- and right-handed surgeons.




101-8173-153D	Curved down	15 cm 6"
101-8173-153U	Curved up	15 cm 6"



STILLE-Ekberg Metzenbaum gracil Diamond SuperCut ergonomisk dissektionssax | STILLE-Ekberg Metzenbaum Delicate Diamond SuperCut Ergonomic Dissecting Scissors

Diamond version for increased durability.




101-8173-153D-C	Curved down	15 cm 6"
101-8173-153U-C	Curved up	15 cm 6"



Metzenbaum STILLE SuperCut XE dissektionssax | Metzenbaum STILLE SuperCut XE Dissecting Scissors

The popular Metzenbaum dissecting scissors with Super Cut edges. Rounded/beveled blades.



101-8175-15	Curved	15 cm 6"
101-8175-18	Curved	18 cm 7"
101-8175-20	Curved	20 cm 8"
101-8175-23	Curved	23 cm 9"
101-8175-26	Curved	26 cm 10.25"
101-8175-30	Curved	30 cm 12"
101-8175-33	Curved	33 cm 13"
101-8176-15	Straight	15 cm 6"
101-8176-18	Straight	18 cm 7"
101-8176-20	Straight	20 cm 8"
101-8176-23	Straight	23 cm 9"
101-8176-26	Straight	26 cm 10.25"



Metzenbaum STILLE Diamond SuperCut dissektionssax | Metzenbaum STILLE Diamond SuperCut Dissecting Scissors

Diamond version of the popular Metzenbaum STILLE SuperCut scissors, for increased durability.



101-8175-15-C	Curved	15 cm 6"
101-8175-18-C	Curved	18 cm 7"
101-8175-20-C	Curved	20 cm 8"
101-8175-23-C	Curved	23 cm 9"
101-8175-26-C	Curved	26 cm 10.25"
101-8175-30-C	Curved	30 cm 12"



Metzenbaum STILLE SuperCut XE ergonomisk dissektionssax | Metzenbaum STILLE SuperCut XE Ergonomic Dissecting Scissors

Metzenbaum Super Cut scissors with ergonomically designed shanks for a more comfortable grip. Rounded/beveled blades.



101-8175-152, 101-8175-182

101-8175-152	Curved	15 cm 6"
101-8175-182	Curved	18 cm 7"



STILLE-Ekberg Metzenbaum SuperCut ergonomisk dissektionssax | STILLE-Ekberg Metzenbaum SuperCut Ergonomic Dissecting Scissors

Ergonomic STILLE-Ekberg version. Available with blades curved up or down, suitable for both left- and right-handed surgeons.



101-8175-153D



101-8175-153U

101-8175-153D	Curved down	15 cm 6"
101-8175-153U	Curved up	15 cm 6"



STILLE-Ekberg Metzenbaum Diamond SuperCut ergonomisk dissektionssax | STILLE-Ekberg Diamond Metzenbaum SuperCut Ergonomic Dissecting Scissors

Diamond version for increased durability.



101-8175-153D-C



101-8175-153U-C


101-8175-153D-C Curved down 15 cm 6"

101-8175-153U-C Curved up 15 cm 6"



Ragnell STILLE dissektions sax | Ragnell STILLE Dissecting Scissors

Dissecting scissors with flat, rounded tips. Often used in plastic surgery. Is also called Killner by some manufacturers. The short models are used for undermining in plastic surgery and the long models are used for further dissection.




101-8378	Straight	12.5 cm 5"
101-8380-12	Curved	12.5 cm 5"
101-8380-15	Curved	15 cm 6"
101-8380-18	Curved	18 cm 7"
101-8380-20	Curved	20 cm 8"
101-8380-23	Curved	23 cm 9"



Ragnell STILLE SuperCut XE dissektions sax | Ragnell STILLE SuperCut XE Dissecting Scissors

Super Cut version of the Ragnell scissors. For precise undermining of flaps.




101-8380-121	Curved	12.5 cm 5"
101-8380-151	Curved	15 cm 6"
101-8380-181	Curved	18 cm 7"
101-8380-201	Curved	20 cm 8"
101-8380-231	Curved	23 cm 9"



Ragnell STILLE Diamond SuperCut dissektionssax | Ragnell STILLE Diamond SuperCut Dissecting Scissors

Diamond version of the Ragnell STILLE SuperCut scissors, for increased durability.




101-8380-121-C	Curved	12.5 cm 5"
101-8380-151-C	Curved	15 cm 6"
101-8380-181-C	Curved	18 cm 7"



STILLE-Ekberg Ragnell SuperCut ergonomisk dissektionssax | STILLE-Ekberg Ragnell SuperCut Ergonomic Dissecting Scissors

Ergonomic STILLE-Ekberg version. Available with blades curved up or down, suitable for both left- and right-handed surgeons.




101-8380-153D	Curved down	15 cm 6"
101-8380-153U	Curved up	15 cm 6"



STILLE-Ekberg Ragnell Diamond SuperCut ergonomisk dissektionssax | STILLE-Ekberg Ragnell Diamond SuperCut Ergonomic Dissecting Scissors

Diamond version for increased durability.



101-8380-153D-C	Curved down	15 cm 6"
101-8380-153U-C	Curved up	15 cm 6"

Ragnell STILLE SuperCut XE ergonomisk dissektionssax | Ragnell STILLE SuperCut XE Ergonomic Dissecting Scissors

Ergo Super Cut versions of Ragnell scissors for precise undermining of flaps and comfortable grip.



101-8380-122




101-8380-152

101-8380-122	Curved	12,5 cm 5"
101-8380-152	Curved	15 cm 6"




Thorek STILLE dissektionssax | Thorek STILLE Dissecting Scissors

Strongly curved dissecting scissors for good visibility. Rounded blades.




101-8470

101-8470	Strongly curved	21 cm 8.25"
----------	-----------------	-------------



STILLE Cooper operations sax | STILLE Cooper Operating Scissors

General purpose utility scissors, often used for sutures, gauze etc. Beveled blades.




101-8042-15	Sharp/blunt	Straight	15 cm 6"
101-8043-15	Blunt/blunt	Straight	15 cm 6"
101-8047-13	Sharp/blunt	Curved	13.5 cm 5.25"
101-8047-15	Sharp/blunt	Curved	15 cm 6"
101-8048-15	Blunt/blunt	Curved	15 cm 6"
101-8048-18	Blunt/blunt	Curved	18 cm 7"




Mayo-Harrington STILLE operations sax | Mayo-Harrington STILLE Operating Scissors

Very useful operating scissors for deep surgery like gynecology and thoracic surgery. Beveled blades.




101-8096-20	Straight	20 cm 8"
101-8096-23	Straight	23 cm 9"
101-8098-20	Curved	20 cm 8"
101-8098-23	Curved	23 cm 9"
101-8098-30	Curved	30 cm 12"



Mayo-Harrington STILLE SuperCut XE operations sax | Mayo-Harrington STILLE SuperCut XE Operating Scissors

Operating Super Cut scissors useful for several surgical disciplines. Beveled blades.



101-8097-20	Straight	20 cm 8"
101-8097-23	Straight	23 cm 9"
101-8099-20	Curved	20 cm 8"
101-8099-23	Curved	23 cm 9"
101-8099-30	Curved	30 cm 12"



Boyd-STILLE tonsillsax | Boyd-STILLE Tonsil Scissors

Scissors with sharp tips and S-shaped shanks for greater visibility. Rounded blades.



101-8420 Curved S-shaped shanks 17 cm 6.75"



Ragnell STILLE operationssax | Ragnell STILLE Operating Scissors

Mammary plastic scissors intended for breast surgery. Rounded blades.



101-8460

101-8460 Angled 17 cm 6.75"



König-STILLE nageluppklippnings sax | König-STILLE Nail-splitting Scissors

Nail-splitting scissors, for cutting up ingrown nails. Beveled blades.



101-8580

101-8580 Straight 15 cm 6"



Stevens-STILLE tenotomisax | Stevens-STILLE Tenotomy Scissors

Tenotomy scissors with rounded, very thin blades and tapered tips. Used in plastic surgery and hand surgery for cutting with extreme precision in small incisions.




101-8166	Straight	12.5 cm 5"
101-8167	Curved	12.5 cm 5"



Stevens-STILLE SuperCut tenotomisax | Stevens-STILLE SuperCut Tenotomy Scissors

Very popular SuperCut XE scissors for precise cutting and dissection of fine tissues. Rounded blades.




101-8166-1	Straight	12.5 cm 5"
101-8167-1	Curved	12.5 cm 5"
101-8167-151	Curved	15 cm 6"
101-8167-181	Curved	18 cm 7"
101-8170-1	Straight	9.5 cm 3.75"
101-8171-1	Curved	9.5 cm 3.75"



Stevens-STILLE Diamond SuperCut tenotomisax
 | Stevens-STILLE Diamond SuperCut Tenotomy Scissors

Diamond version for increased durability.




101-8167-1-C	Curved	12.5 cm 5"
101-8167-151-C	Curved	15 cm 6"
101-8167-181-C	Curved	18 cm 7"



Stevens-STILLE SuperCut XE ergonomisk tenotomisax | Stevens-STILLE SuperCut XE Ergonomic Tenotomy Scissors

Ergonomic design of the popular Stevens SuperCut XE tenotomy scissors.




101-8166-2	Straight	12.5 cm 5"
101-8167-2	Curved	12.5 cm 5"
101-8167-152	Curved	15 cm 6"
101-8167-182	Curved	18 cm 7"



Stevens-STILLE Diamond SuperCut ergonomisk tenotomisax | Stevens-STILLE Diamond Ergonomic SuperCut Tenotomy Scissors

Diamond version for increased durability.



101-8167-2-C

101-8167-2-C Curved 12.5 cm 5"



STILLE-Ekberg Stevens SuperCut ergonomisk tenotomisax | STILLE-Ekberg Stevens Ergonomic SuperCut Tenotomy Scissors

Ergonomic STILLE-Ekberg version. Available with blades curved up or down, suitable for both left- and right-handed surgeons.



101-8167-153D



101-8167-153U

101-8167-153D Curved down 15 cm 6"

101-8167-153U Curved up 15 cm 6"



STILLE-Ekberg Stevens Diamond SuperCut ergonomisk tenotomisax | STILLE-Ekberg Stevens Diamond Ergonomic SuperCut Tenotomy Scissors

Diamond version for increased durability.



101-8167-153D-C



101-8167-153U-C


101-8167-153D-C Curved down 15 cm 6"

101-8167-153U-C Curved up 15 cm 6"




STILLE SuperCut XE kärlsax | STILLE SuperCut XE Vascular Scissors

Super Cut scissors with angled blades for better reach of vessels etc.
Used in vascular surgery.




101-8200-25	Angled 25°	18 cm 7"
101-8200-45	Angled 45°	18 cm 7"
101-8200-60	Angled 60°	18 cm 7"
101-8200-90	Angled 90°	17 cm 6.75"
101-8200-125	Angled 125°	17 cm 6.75"



STILLE SuperCut XE kärlsax | STILLE SuperCut XE Vascular Scissors

Extra delicate version of the popular STILLE Vascular Scissor with shorter and thinner blades.



101-8201-25	Angled 25°	18 cm 7"
101-8201-45	Angled 45°	18 cm 7"
101-8201-60	Angled 60°	18 cm 7"
101-8201-90	Angled 90°	17 cm 6.75"
101-8201-125	Angled 125°	17 cm 6.75"



STILLE-Lund kärlsax | STILLE-Lund Vascular Scissors


Vascular scissors with one sharp and one probe-pointed blade for precise handling of fine blood vessels. Popular model in pathology procedures.

101-8490	1 probe-pointed blade	Straight	15 cm 6"
101-8490-2	1 probe-pointed blade	Ergo	Straight 15 cm 6"



**Fomon STILLE nässax |
Fomon STILLE Nasal Scissors**

Scissors for rhinoplasty. Double-beveled blades.




101-8025

101-8025	Curved	12,5 cm 5"
----------	--------	------------

**Joseph STILLE nässax |
Joseph STILLE Nasal Scissors**

Nasal scissors for the lower lateral cartilage. Useful in dissecting and cutting tough tissue. Beveled blades.



101-8110



101-8115

101-8110	Straight	15 cm 6"
101-8115	Curved	15 cm 6"

**Joseph STILLE SuperCut XE nässax | Joseph
STILLE SuperCut XE Nasal Scissors**

Super Cut nasal scissors for the lower lateral cartilage. Useful in dissecting and cutting tough tissue. Beveled blades.



101-8110




101-8115

101-8110-1	Straight	15 cm 6"
101-8115-1	Curved	15 cm 6"

Peck-Joseph STILLE nässax |
Peck-Joseph STILLE Nasal Scissors


Double-beveled scissors for rhinoplasty.



101-8120	Sharp	Straight	15 cm 6"
101-8122	Blunt	Straight	15 cm 6"
101-8123	Sharp	Curved	15 cm 6"
101-8125	Blunt	Curved	15 cm 6"

Peck-Joseph STILLE SuperCut XE nässax |
Peck-Joseph STILLE SuperCut XE Nasal Scissors

Double-beveled scissors for rhinoplasty.



101-8120-1	Sharp	Straight	15 cm 6"
101-8123-1	Sharp	Curved	15 cm 6"

STILLE SuperCut XE nässax | STILLE SuperCut XE Nasal Scissors

SuperCut XE nasal scissors for excision of cartilage.



101-8411

101-8411 Straight Angled shanks 18 cm 7"



STILLE Fomon SuperCut XE nässax | STILLE Fomon SuperCut XE Nasal Scissors

Serrated SuperCut XE scissors for cutting dorsal cartilage. Angled shanks for better visibility.




101-8425

101-8425 Serrated blades Angled shanks 14 cm 5.5"



STILLE Cottle SuperCut XE nässax | STILLE Cottle SuperCut XE Nasal Scissors

SuperCut XE scissors with heavy blades for cutting dorsal cartilage. Angled shanks for better visibility. Serrated blades.





101-8604

101-8604 Serrated blades Angled shanks 16 cm 6.25"



Gorney STILLE SuperCut XE Face-lift sax | Gorney STILLE SuperCut XE Face-lift scissors



Super Cut double-beveled face-lift scissors with flat tips and semi-sharp outer blades. Two models with serrated blades to prevent tissue slippage.



101-8430-171	Straight		17 cm 6.75"
101-8431-171	Curved		17 cm 6.75"
101-8433-20	Curved		20 cm 8"
101-8435-20	Curved	Serrated	20 cm 8"
101-8436-20	Straight	Serrated	20 cm 8"

Gorney STILLE SuperCut XE ergonomisk Face-lift sax | Gorney STILLE SuperCut XE Ergonomic Face-lift scissors

Ergonomic design to avoid fatigue in the hand. One model with serrated blades to prevent tissue slippage. Flat tips and semi-sharp outer blades.



101-8430-172	Straight		17 cm 6.75"
101-8431-172	Curved		17 cm 6.75"
101-8433-202	Curved		20 cm 8"
101-8435-202	Curved	Serrated	20 cm 8"

Aston STILLE SuperCut XE Face-lift sax | Aston STILLE SuperCut XE Face-lift scissors


Double-beveled SuperCut XE face-lift scissors with semi-sharp outer blades for undermining flaps.



101-8440-201



101-8441-201



101-8442-171

101-8440-201	Straight		20 cm 8 in
101-8441-201	Curved		20 cm 8 in
101-8442-171	Curved	Serrated	17 cm 6.75 in



Davis STILLE SuperCut XE Face-lift sax | Davis STILLE SuperCut XE Face-lift Scissors

SuperCut XE face-lift scissors with tapered blades and sharp tips for precise cutting and dissection.



101-8452

101-8452	Curved	18 cm 7"
----------	--------	----------



**Rees STILLE SuperCut XE Face-lift sax |
Rees STILLE SuperCut XE Face-lift scissors**

SuperCut XE face-lift scissors, double-beveled with semi-sharp outer blades for undermining skin flaps.




101-8456-171

101-8456-171	Curved	17 cm 6.75"
--------------	--------	-------------



**Matarasso STILLE SuperCut XE Face-lift sax |
Matarasso STILLE SuperCut XE Face-lift scissors**

SuperCut XE face-lift scissors with double-beveled semi-sharp outer blades.



101-8458

101-8459-20

101-8458-17	Curved	17 cm 6.75"
-------------	--------	-------------

101-8458-20	Curved	20 cm 8"
-------------	--------	----------

101-8459-20	Straight	20 cm 8"
-------------	----------	----------



**Matarasso STILLE SuperCut XE ergonomisk
Face-lift sax | Matarasso STILLE SuperCut XE
Ergonomic Face-lift scissors**

Ergonomic version of the Matarasso SuperCut XE face-lift scissors with double-beveled semi-sharp outer blades.



101-8458

101-8458-172	Curved	17 cm 6.75"
--------------	--------	-------------

101-8458-202	Curved	20 cm 8"
--------------	--------	----------

Littler STILLE trådföringssax |
Littler STILLE Suture carrying scissors


Suture carrying scissors. Rounded blades.



101-8164

101-8164 Curved 12,5 cm 5"






STILLE sutursax | STILLE Suture scissors

Suture removal scissors.



101-8584

101-8584 Curved 13.5 cm 5.25"



Spencer STILLE sutursax | Spencer STILLE Suture scissors

Popular suture scissors for many procedures. Delicate model.

101-8586-9 9 cm 3.5"



Lister STILLE TC bandagesax | Lister STILLE TC Bandage Scissors

Bandage scissors with TC-inserts for improved durability.


101-8595-18 TC-inserts 18 cm 7"



**Littauer STILLE sutursax |
Littauer STILLE Suture scissors**

A somewhat larger model than Spencer. Used in many procedures.

101-8597-14 Delicate pattern 14 cm 5.5"



**Universal STILLE ståltrådssax |
Universal STILLE Wire Cutting scissors**

Wire cutting scissors with TC-inserts. For soft wire up to 1.2 mm,
for hard wire up to 0.8 mm.

101-8815-12 TC 12 cm 4.75"